

Product datasheet for MR218630

Negr1 (NM_001039094) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Negr1 (NM_001039094) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Negr1
Synonyms:	5330422G01Rik; KILON; Ntra
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR218630 representing NM_001039094 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGTGCTCCTGGCGCAGGGCGCTTGTGCTCGAACCAGTGGCTGGCGGGGTGCTTCTGAGCCTGTGCT
CTTGCCTCCCGCCGGGAGAGCGTGGACTTCCCTTGGCGGCCGTGGACAACATGCTGGTGAGGAAAGG
TGACACAGCGGTGCTCAGGTGTTACTTGAAGATGGAGCATCAAAGGGTGCCTGGCTCAACAGGTCAAGT
ATCATTGCTGGAGGTGACAAGTGGTCAAGTGGACCCTCGAGTTCCATTTCCACATTGAATAAAAGAG
ACTACAGCCTCCAGATACAGAATGTTGATGTAACAGATGATGGCCCGTACACCTGTTCTGTGCAGACCCA
ACACACACCACGGACAATGCAGGTTCACTCACTGTGCAAGTCCCACCGAAAATATACGACATCTCAAAT
GACATGACCATCAATGAAGGAACCAATGTACCCTTACTTGTGGCCACTGGGAAGCCAGAGCCTGTCA
TTTCTGGAGGCATATCTCCCATCAGCAAAACCATTTGAAAATGGACAATATTTGGACATTTATGGAAT
TACAAGGGACCAGGCTGGGGAGTATGAATGCAGTGCAGAAAACGATGTGTGTTCCAGACGTGAAGAAA
GTGAGAGTGATCGTGAACCTTTCGCCTACAATTCAGGAAATTAATCTGGCACTGTGACCCTGGACGCA
GTGGACTGATACGATGCGAGGGTGCAGGTGTGCCGCCCGCCCTTCGAGTGGTACAAAGGAGAGAAGAG
GCTCTTCAATGGCCAACAAGGAATTATCATTGAGAAATTTAGCACAAGATCCATCCTCACAGTGACCAAT
GTGACGCAGGAGCACTTCGGCAACTATACTTGTGTGGCTGCCAACAAGTTGGGTACAACCAACGCGAGCC
TGCCCTCAACCTCCAAGTACAGCCAGTATGGAATTAAGTGGAGTGCCTGCGACCTCTTCTCCTGCTG
GAGCCTTGCCTGACACTATCTTCTGTCATCAGTATATTCTACCTGAAGAATGCTATCCTCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR218630 representing NM_001039094
 Red=Cloning site Green=Tags(s)

MVLLAQQGACCSNQWLAAVLLSLCCLPAGQSVDFPWAAVDNMLVRKGD AVLRCYLEDGASKGAWLNRSS
 IIFAGGDKWSVDPRVISISTLNKRDSLQIQNVDVTDDGPYTC SVQTQHTPRMQVHLTVQVPPKIYDISN
 DMTINEGTNVTLTCLATGKPEPVISWRHISPSAKPFENGQYLDIYGITRDQAGEYEC SAENDVSFPDVKK
 VRVIVNFAPTIQEIKSGVTTPGRSGLIRCEGAGVPPPAFEWYKGEKRLFNGQQGIIIQNFSTRSILTVTN
 VTQEHFGNYTCVAANKLGTNASLPLNPPSTAQYGITGSACDLFSCWSLALTLSSVISIFYLKNAILQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001039094

ORF Size: 1044 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039094.3](#)

RefSeq Size: 4985 bp

RefSeq ORF: 1047 bp

Locus ID: 320840

UniProt ID: [Q80Z24](#)

Cytogenetics: 3 H4

MW: 38.4 kDa

Gene Summary: May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain.[UniProtKB/Swiss-Prot Function]