

## Product datasheet for **MR218651L3V**

### Nt5c2 (NM\_001164363) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nt5c2 (NM_001164363) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nt5c2
Synonyms:	2010002I23Rik; cN-II; CnII; ENSMUSG00000025039; Gm9751; Gmp; Nt5b; Pnt5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164363
ORF Size:	1758 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218651).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164363.1</a>
RefSeq Size:	4181 bp
RefSeq ORF:	1761 bp
Locus ID:	76952
UniProt ID:	<a href="#">Q3V1L4</a>
Cytogenetics:	19 C3



[View online »](#)

**Gene Summary:**

May have a critical role in the maintenance of a constant composition of intracellular purine/pyrimidine nucleotides in cooperation with other nucleotidases. Preferentially hydrolyzes inosine 5'-monophosphate (IMP) and other purine nucleotides (By similarity). [UniProtKB/Swiss-Prot Function]