

Product datasheet for MR218704

Dmrtc1a (NM_027591) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dmrtc1a (NM_027591) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dmrtc1a
Synonyms: 1700030H14Rik; 4921520P21Rik; Dmrt8.1; DMRTC1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR218704 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGCGGTCTTCGGGTTCCCGGGAACCTCGTAAGGACTTTTCTCCTGTGACCACCTCCAAGAAGGAGG
 AGGCTACCCCATTAAGAGACGCCTGGTTGAGAGACATAAAAGGACTATGGCTGCTGCTCATCCATCACA
 CATGCACGTCAAGAAATGGCTGTAGAAGAAGGAGTTCGCACTGGAAAGAACAACACTGTGCAGCAGATCCAG
 GCCCAGGTTGACACAGCAACCCAGGAGGAGAGCTCCCAAGGCCCTGTGCTGCTGAACCAGCACCCAGAAA
 CAACATCTGTGCCATATACTCCAGAGACTGTGGGCAACAACATGATGGTGTCCCTTCCTGGGAGCCCCA
 TGGGACCTCTGCCATGCCAAGCATGTGTCCAAGCCTGATTCCTCAGCCCTGTGCCACCACTGACCCATATG
 CTAAGTCAACACAGGTCATGGGCCCCAGTGCCTCCAACCAGGCTTCAGTTTCTGCTACCTTAGAGTGGC
 AGGAGATGCTGGAGGCTGCTGAAGCTCTGCTGGCTCTGAAAAATTCGTCTCAGACACGTCACCAGCCCTG
 TGGCATGCCAGTGAAGAGAGAGCTAAGCGGAGCATCACGTATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218704 protein sequence
 Red=Cloning site Green=Tags(s)

MQRSSGSREPRKDFSPVTTSKKEEATPLKRRLVERHKRTMAAAHPSMHVKKLAVEEGVRTGKNTVQQIQ
 AQVDTATQEESQGPVLLNQHPETTSVPYTPETVGQQLMVSLPGEPHGTSAMPSMCPSLILQPCATDPM
 LLQPQVMGPSASNQASVSATLEWQEMLEAAEALLALKNSSQTRHQPCGMPVKRELSGASRM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_027591

ORF Size: 606 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_027591.3](#), [NP_081867.1](#)

RefSeq Size: 1247 bp


RefSeq ORF: 606 bp

Locus ID: 70887

UniProt ID: [Q9D9R7](#)

Cytogenetics: X D
 MW: 21.9 kDa

Product images:



Circular map for MR218704