

Product datasheet for MR218717

Nkx2-9 (NM_008701) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nkx2-9 (NM_008701) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nkx2-9
Synonyms:	Nkx-2.9; Nkx2-8; Nkx2.9; tinman
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR218717 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGGCCACCTCTGGACGCTCGGCTTCACGGTGGCAGCCTCCTGAATTTACCTGAACAGGATGCGAAGC
CTAGGGTGAGGCGGAGCAACAGACGTGCGTTCCACAGACCGCTGCCTGGCTGGAGTCGGAGTGCAGCCA
CTACCTGTCTCGGACGAGAGCGGCTGGAGACAAGCCCTGCAGACTCGTCGAGCTGGCTTCGCTCAGG
CGGGAGTCCCCGGCTCCGACCCGGAGAAGAGGAGGAAGCGGCGGGTGTGTTCTCAAGGCCAGACCT
TGGAGTTGGAGCGGCGCTTTCGGCAGCAGCGCTACCTGTGGCGCCGGAGCGGGAACAGCTGGCTCGACT
CCTGCGCCTCACGCCACGCAGGTCAAGATCTGGTTCAGAACCCCGCTACAAGCTGAAGCGAGGGCGA
GCGCCTGGGATCACAGAGCCTTCGGATATGGCAGCATCATCTGATCTGCACGCCGCTCCCGGCTGCTGC
GCCGCTAGTGGTACCGTGTGGTTCACGACCGGCCACCGAGCAACAACGGCAGAGGTGAGGGGACCTC
TGCGGTGCCCCAGGACAAGTGCAGCGCACGCTGGCCACTGCGTGTCCGGTGCCAGGCTACACCGCCTTC
GGACCGGCTCGGCGCTTGGCCTTCCAGCCTACCAGCACTTAGCACCACCAGCTCTGGTCTCCTGGA
ACTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR218717 protein sequence
Red=Cloning site Green=Tags(s)

MATSGRLGFTVRSLLNLPEQDAKPRVRREQQTCVPQTAAWLESECSHYLSSDESGLTSPADSSQLASLR
 RESPGSDPEKRRRRLVLF SKAQ TLELERRFRQRYLSAPEREQLARLLRLTPTQVKIWFQNHRYKLRGR
 APGITEPSDMAASDLHAAPGLLRVVPVLVHDRPPSNNGREGTSAVPQDKCSARLATACPVPGYTAF
 GPGSALGLFPAYQHLAPPALVSNW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_008701

ORF Size: 708 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008701.1](#), [NM_008701.2](#), [NP_032727.2](#)

RefSeq Size: 1260 bp

RefSeq ORF: 708 bp

Locus ID: 18094


UniProt ID: [O70584](#)

Cytogenetics: 12 C1

MW: 26.1 kDa

Gene Summary: Possible role in the specification of a distinct subset of neurons.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR218717