

Product datasheet for MR218742

Cnih3 (NM_001160211) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cnih3 (NM_001160211) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cnih3
Synonyms: 2900075G08Rik; CNIH-3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR218742 representing NM_001160211
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCCTTCACTTTTCGCTGCTTTCTGCTACATGCTGTCCCTGGTGTCTGCGCTGCGCTCATCTTTTTCG
CTATCTGGCACATAATCGCCTTTGACGAGCTAAGAACGGATTTCAAGAGTCCCATAGACCAGTCAACCC
GGTTCATGCGAGAGAACGGCTGAGGAACATCGAACGCATCTGCTTCTTCTGAGAAAGCTGGTGCTGCC
GAGTACTCTATCCACAGCCTCTTCTGCATCATGTTCTGTGCGCTCAGGAGTGGCTCACCCCTGGGACTCA
ATGTTCTCTGCTTTTCTATCACTTCTGGAGGTATTTCCACTGTCCAGCAGATAGCTCAGAGCTAGCCTA
CGACCCACCAGTTGCATGAATGCGGACACCTTGAGCTACTGCCAGAAGGAGGCCTGGTGTAAAGCTGGCC
TTCTACCTCTCTCTTCTTCTACTACCTGTACTGCATGATTTACACCTTAGTGAGCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218742 representing NM_001160211
Red=Cloning site Green=Tags(s)
MAFTFAAFCYMLSLVLCALIFFAIWHIIAFDELRTDFKSPIDQCNPVHARERLRNIERICFLLRKLVLV
EYSIHSFLFCIMFLCAQEWLTLGLNVPLLFYHFWRYFHCPADSSSELAYDPPVVMNADTLSYCQKEAWCKLA
FYLLSFFYYLYCMIYTLVSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001160211

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001160211.1</u> , <u>NP_001153683.1</u>
RefSeq Size:	2670 bp
RefSeq ORF:	483 bp
Locus ID:	72978
UniProt ID:	<u>Q6ZWS4</u>
Cytogenetics:	1 H4
MW:	19.4 kDa
Gene Summary:	Regulates the trafficking and gating properties of AMPA-selective glutamate receptors (AMPA receptors). Promotes their targeting to the cell membrane and synapses and modulates their gating properties by regulating their rates of activation, deactivation and desensitization (By similarity).[UniProtKB/Swiss-Prot Function]