

Product datasheet for **MR218768L4V**

Atp7b (NM_007511) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp7b (NM_007511) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atp7b
Synonyms:	Atp7a; tx; WND
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007511
ORF Size:	4386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218768).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007511.2 , NP_031537.2
RefSeq Size:	4711 bp
RefSeq ORF:	4389 bp
Locus ID:	11979
UniProt ID:	Q64446
Cytogenetics:	8 10.78 cM
Gene Summary:	Copper ion transmembrane transporter involved in the export of copper out of the cells, such as the efflux of hepatic copper into the bile.[UniProtKB/Swiss-Prot Function]



[View online »](#)