

Product datasheet for **MR218774L4V**

Cd180 (NM_008533) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd180 (NM_008533) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd180
Synonyms:	F630107B15; Ly-78; Ly78; RP105
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008533
ORF Size:	1983 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218774).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008533.2 , NP_032559.2
RefSeq Size:	2326 bp
RefSeq ORF:	1986 bp
Locus ID:	17079
UniProt ID:	Q62192
Cytogenetics:	13 D1
Gene Summary:	May cooperate with MD-1 and TLR4 to mediate the innate immune response to bacterial lipopolysaccharide (LPS) in B-cells. Leads to NF-kappa-B activation. Also involved in the life/death decision of B-cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)