

## Product datasheet for MR218805

### Bst1 (NM\_009763) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Bst1 (NM\_009763) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Bst1  
**Synonyms:** 114/A10; A530073F09; BP-3; Bp3; Bsta1; CD157; Ly65  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR218805 representing NM\_009763  
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTGTCCAGGGAGGCCTGTTGTCTCTGTGGCTGTGGCTATGGCTGTCCTTGCTGACAGTGCTTCTGG  
 GGGCAAGAGCGCGGTGGACAGGGGAGGGTACCCTCCTCACCTCCAGAGCATCTTCTGGGCCGCTGCGC  
 AGAGTACACCACGCTGCTGAGCCTGGGAACAAAACTGCACAGCCATCTGGGAGCCTTCAAGGGGGTG  
 CTGGACAAGGACCCTGCTCTGTGCTCCCCTCAGACTATGATCTCTTCAACCTCTCCAGGCACCCCA  
 TTCTAGGACAAGTCACTGTTCTGGGAGAATAACCACCTGCTCGTTATGAGCTATGGGAGAACACCCG  
 CCGCCTCGTGGCCCTGTGCGACGTGCTGTATGGCAAAGTTGGCGATTTCTTGAGCTGGTGTGACAGGAA  
 AATGCCTCTGGACTTGACTACCAGTCTGCCCCACATCAGAAGACTGTGAAAACAACGCCGTGGACTCCT  
 ATTGAAAAGTGCATCCATGCAGTATCAAGAGACAGTTCAGGGGTGATCAACGTGATGCTGAATGGCTC  
 GGAGCCAAAAGGAGCCTATCCCACGAGAGGGTTTTTTCAGATTTTGAATCCCCTACTTACAAAAGGAT  
 AAAGTACAAGAATAGAGATCTGGGTTATGCATGACGTTGGAGGGCCCAATGTGGAATCCTGTGGCGAAG  
 GCAGTGTAAAACTCTGGAAGACAGATTGGAGGCCCTGGGTTCCAGCACAGCTGTATCAATGACTACCG  
 ACCAGTGAAGTTCTTAATGTGTGTGGACCACAGCACTACCCCGACTGTATCATGAACTCAGCCTCAGCA  
 TCTATGCGGAGGGAGAGCGCATCTTTCATGCAATAGGAGACGCCAGCCTTCTCATTCTCTTAGTAG  
 CTTTGGCTTCAAGTTCTCAGGCG


**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**



**ACCN:** NM\_009763

**ORF Size:** 933 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_009763.3](#), [NP\\_033893.2](#)

**RefSeq Size:** 2336 bp

**RefSeq ORF:** 936 bp

**Locus ID:** 12182


**UniProt ID:** [Q64277](#)

**Cytogenetics:** 5 23.84 cM

**MW:** 35 kDa

**Gene Summary:** Synthesizes the second messengers cyclic ADP-ribose and nicotinate-adenine dinucleotide phosphate, the former a second messenger that elicits calcium release from intracellular stores. May be involved in pre-B-cell growth.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR218805