

Product datasheet for **MR218807L3V**

A4galt (NM_001170954) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	A4galt (NM_001170954) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	A4galt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001170954
ORF Size:	1080 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218807).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001170954.1 , NP_001164425.1
RefSeq Size:	2145 bp
RefSeq ORF:	1047 bp
Locus ID:	239559
UniProt ID:	Q67BJ4
Cytogenetics:	15 E1



[View online »](#)

Gene Summary:

The protein encoded by this gene catalyzes the transfer of galactose to lactosylceramide to form globotriaosylceramide, which has been identified as the P(k) antigen of the P blood group system. The encoded protein, which is a type II membrane protein found in the Golgi, is also required for the synthesis of the bacterial verotoxins receptor. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2010]