

Product datasheet for **MR218835**

Nxt2 (NM_001161430) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nxt2 (NM_001161430) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nxt2
Synonyms: 6330587F24Rik; AI787442; P15-2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR218835 representing NM_001161430
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCTTGGCAGTGAATTTTAAAACCTTATGTGGATCAGGCATGTAGAGCTGCCGAGGAATTTGTCAATA
TTTACTATGAGACAATGGACAAAAGACGACATGCACCTGTCAGGCTGATTTGGACAAAGCCACTCTAAT
CTGGAATGGAAATGTTGTTACAGGGCTGGAAGCCCTTGCTAATTTTTTGAGATGTTGCCTCCAGCGAA
TTCCAGATCAATATGTTAGATTGCCAACCGTTCATGAGCAAGCTACCCAGTGCCAAACTACAGTGCTCG
TTGTGACCAGTGGAGTTGTGAAGTTTGTGAAACAAGCAGCACTTCTTTAACCAGAACTCCTGCTGAC
TGCGCAGAGCACTCCTAACAGCACTGTGTGGAAGATTGCAAGTGATTGCTCCGTTTTCAAGATTGGCT
AGTATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218835 representing NM_001161430
Red=Cloning site Green=Tags(s)

MALAVNFKTYVDQACRAAEFVNIYYETMDKRRHALVRLYLDKATLIWNGNVVTGLEALANFFEMLPSSSE
FQINMLDCQPVHEQATQCQTTLVVTSGVVKFDGNKQHFNFQNFLLTAQSTPNSTVWKIASDCFRFDWA
SI

TRTRPLEQK**L**ISEEDLAANDILDYKDDDDK**V**


Chromatograms: https://cdn.origene.com/chromatograms/mm9053_h07.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001161430

ORF Size: 426 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161430.2](#), [NP_001154902.1](#)

RefSeq Size: 2455 bp

RefSeq ORF: 429 bp

Locus ID: 237082


UniProt ID: Q3UNA4

Cytogenetics: X F2

MW: 16.3 kDa

Gene Summary: Regulator of protein export for NES-containing proteins. Also plays a role in mRNA nuclear export (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR218835