

Product datasheet for MR218870L3

Idh1 (NM_001111320) Mouse Tagged Lenti ORF Clone


Product data:

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

Product Type:	Expression Plasmids		
Product Name:	ldh1 (NM_001111320) Mouse Tagged Lenti ORF Clone		
Tag:	Myc-DDK		
Symbol:	ldh1		
Synonyms:	Al314845; Al788952; E030024J03Rik; ld-1; ldh-1; ldpc		
Mammalian Cell Selection:	Puromycin		
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)		
E. coli Selection:	Chloramphenicol (34 ug/mL)		
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218870).		
Restriction Sites:	Sgfl-Mlul		
Cloning Scheme:	Cloning sites used for ORF Shuttling:		
	GCG ATC GCC ATG // NNN ACG CGT		




* The last codon before the Stop codon of the ORF.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Plasmid Map:



ACCN:	NM_001111320	
ORF Size:	1242 bp	
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>	
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.	
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).	
Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. 	
RefSeq:	<u>NM 001111320.1, NP 001104790.1</u>	
RefSeq Size:	2292 bp	
RefSeq ORF:	1245 bp	

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	ldh1 (NM_001111320) Mouse Tagged Lenti ORF Clone – MR218870L3	
Locus ID:	15926	
UniProt ID:	<u>088844</u>	
Cytogenetics:	1 32.91 cM	

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US