

Product datasheet for MR218932L3

Uts2 (NM_011910) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Uts2 (NM_011910) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Uts2
Synonyms:	pre; prepro-Ull; Ucn; Ucn2; Ull
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218932).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_011910
ORF Size:	372 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011910.2 , NP_036040.1
RefSeq Size:	538 bp
RefSeq ORF:	372 bp
Locus ID:	24111
UniProt ID:	Q9QZQ3
Cytogenetics:	4 E2
Gene Summary:	This gene encodes a member of the urotensin family of peptide hormones. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional hormone before secretion into the plasma. Mice lacking the encoded protein have a significantly decreased low density lipoprotein cholesterol profile and hepatic steatosis that is consistent with decreased hepatocyte de novo cholesterol synthesis and apolipoprotein B secretion. [provided by RefSeq, Sep 2016]

Product images:



Circular map for MR218932L3