

## Product datasheet for **MR218939L4V**

### **Pck2 (NM\_028994) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pck2 (NM_028994) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pck2
Synonyms:	1810010O14Rik; 9130022B02Rik; PEPCK-M
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028994
ORF Size:	2001 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218939).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028994.2</a> , <a href="#">NP_083270.1</a>
RefSeq Size:	3400 bp
RefSeq ORF:	1923 bp
Locus ID:	74551
UniProt ID:	<a href="#">Q8BH04</a>
Cytogenetics:	14 28.19 cM
Gene Summary:	Catalyzes the conversion of oxaloacetate (OAA) to phosphoenolpyruvate (PEP), the rate-limiting step in the metabolic pathway that produces glucose from lactate and other precursors derived from the citric acid cycle.[UniProtKB/Swiss-Prot Function]



[View online »](#)