

Product datasheet for **MR218948**

Odf2l (NM_025714) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Odf2l (NM_025714) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Odf2l
Synonyms:	4733401D09Rik; 9630045K08Rik; D3Ert250e
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR218948 representing NM_025714
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGAGATGCCTACTAGTGATGGAAGTCTTTCAGAGGCGCTGATTCTTATCTCATCCGGGAGAGAGGAT
 CTGCTCCACAGTGTGCCCTGAGAAGCACCTCAGCAGGCTCAAGCAAGACGTCGCAAAATACAAAGATGGA
 GTTGAAAGCAACGCTTAAGGAAGCTCAGTTGGCATCCTGTTCAAGTTGAATTACTTTTGCCATTATTTAAG
 AACACTGTGGAAGGAATTAGTCTAGAAAACGCTAATCTTCTGCATCCAGTTTGAAGAAGATTTTGAAC
 AGAACGACATATTATCAAAAGAATTAGACACTTTTAACCGTGTGAAATTAGCCTTAGAACATCTTATTA
 ACAGACTGACTATGAACAAATAGGAGACAGCTTGCTATGCTTGTGAAGGATCTCAGTGACAATGAAAGT
 GAGAACAGAAACCTTGAGGAAAAAGTACTTGAAGGAGACCTATATCCGTGAACCTCTCGTGTCTGTTTC
 AGAATGAAAGGAAAGTCTTTGAAAGCAAACAGGCTTTCTCAATCAGTCAAAGTAGTCACGATCGTCT
 ACAGCGCCAAATTCAGAAGAGGGAAGCTGAGAATGAAAGCTAAAAGAACATGTACAGAGCTTGAAACC
 CAAATAGCCAAGTGGAACCTTGCAAGTGAAAATGAATAAGCAAGAGGCTGTAGCCATCAAGGAAGCAAGTC
 GGCAAAAGGCTGTAGCTCTGAAAAAGCATCGAAAGTGTACAGACAAAGGCTGCGACACTTCACAGGGGA
 TATCGAGCGCCTCGCCTCCCAAGTCAGAGATCAGGAAGCCAAGCTGTCTGAAACAGTTTCAGCTTCCAGT
 GACTGGAAAGTCAGTTCGAGAAATTTGCCATAGAAAAAAGTCTGAGCTGGAAGTGCAGATCGAAACATGA
 AAAAGCAAATCGCTAACCTTTTGAAGACTTGAGAAAAATGGAAGTCTGGAAGAACTCATGTGAAGA
 AATTCCTAGAAAACCTCACTCACTGGAAGATGAAATGAAGCCCTAAATATTGAGAATGTGAAATTAAG
 GGTACACTTGATGCCTTGAAGGATGAAGTTGCTTCTGTTGAAATGAAGTGTAGAACTGCAAGAAGTAG
 AAAACGACAAAAACCTCGTCAAGGATACAGGACTCAGGTACAAAAGTTGCAAGAAGCAGCTGAAAT
 GGTGAAAAGTAGATGCAAAAATTTACTACATGAAAATAACCTAATAATAACGAATAAAAAACAAAAGTTA
 GAGAAGGTTGATGGAACCAAAGCCTTCTGACCAAGCTTCTCTTGAAGGAAAGAAATCATCTCATTCACT
 TAAAGTGTGAGAACCTTAAAGAGAACTAGAGCAGATGGATGCAGAGAATAAAGAGCTTGAGAAGAAGCT
 GGCCGACCAAGAGGAGTGTCTCAAGCACAGTGACCTGGAGCTGAAAGAGAAGGCTGCAGAGTACACAGCC
 CTGTCCAGGCAACTGGAAGCCGCTTTGGAAGAAGGAAGACAAAAGGTTTCAGAAGAAGTAGAGAAAATGT
 CATCTAGAGAGAGGCTTTACAAATTAATAATTAGATCTGGAAGCTGAGCTTAGAAAGAAAAATGAAGA
 ACAAACCAACTTGTGACAAGATGAACACTAAACCCAGCACCAGGCTATCTGCTTGAAGGAAATACAG
 CACAGCCTTGAGAAGTCGGAGACTCGGAATGAGAGCATTAGAATTACTGCAGTTTCTTCAATTTCTT
 ATGTAACGATGTTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR218948 representing NM_025714
 Red=Cloning site Green=Tags(s)


MEMPTSDGSLSEGLISYLIRERGSAPQCAPEKHL SRLKQDVANTKMEKATLKEAQLASCSVELLLPLFK
 NTVEGISLENANLSASSLKKIFEQNDILSKELDTFNVRKLAL EHLIKQTDYEQIGDSSLCLLKDLSDNES
 ENRNLEEKVLEKETYIRELSCLFQNEKESALKANRLSQSVKVVHDLRQRQIKREAENEKLEHVQSLET
 QIAKWNLVQVMNKQEAVAIKEASRQKAVALKKASKVYRQLRHFTGDIERLASQVRDQEAKLSETVSASS
 DWKSQFEKIAIEKTELEVQIETMKKQIANLLEDLRKMETHGKNSCEEILRKLHSL EDENEALNIENVKLL
 GTLDALKDEVASVENELVELQEVEKRQKTLVEGYRTQVQKLQEAEMVKSRCNLLHNNLIITNKNKKL
 EKVDGNQSLTKLSLEEENHLIQLKCENLKEKLEQMDAENKELEKKLADQEECLKHSDELKEKAAEYTA
 LSRQLEAALEEGRQKVSEEEVKMSSRERLQIKILDLEAELRKKNEEQNLVDKMNTKTQHQAICLKEIQ
 HSLEKSETRNESIKNYLQFLQISVYTMFR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM 025714

ORF Size: 1767 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_025714.4](#), [NP_079990.2](#)

RefSeq Size: 2291 bp

RefSeq ORF: 1770 bp

Locus ID: 52184


UniProt ID: [Q9D478](#)

Cytogenetics: 3 69.36 cM

MW: 68.5 kDa

Gene Summary: Acts as a suppressor of ciliogenesis, specifically, the initiation of ciliogenesis.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR218948