

Product datasheet for **MR218961L3V**

Tktl1 (NM_031379) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tktl1 (NM_031379) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tktl1
Synonyms:	C79967
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031379
ORF Size:	1785 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218961).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031379.2 , NP_113556.2
RefSeq Size:	2485 bp
RefSeq ORF:	1788 bp
Locus ID:	83553
UniProt ID:	Q99MX0
Cytogenetics:	X A7.3
Gene Summary:	Catalyzes the transfer of a two-carbon ketol group from a ketose donor to an aldose acceptor, via a covalent intermediate with the cofactor thiamine pyrophosphate. [UniProtKB/Swiss-Prot Function]



[View online »](#)