

Product datasheet for **MR218967L3V**

Dfnb59 (Pjvk) (NM_001080711) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Pjvk
Synonyms:	Dfnb59; Gm1001; Pjvk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001080711
ORF Size:	1056 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR218967).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001080711.2 , NP_001074180.1
RefSeq Size:	1217 bp
RefSeq ORF:	1059 bp
Locus ID:	381375
UniProt ID:	Q0ZLH2
Cytogenetics:	2 C3



Gene Summary: Essential in the activity of auditory pathway neurons.[UniProtKB/Swiss-Prot Function]