

Product datasheet for **MR218973L3V**

Rgl3 (NM_023622) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rgl3 (NM_023622) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rgl3
Synonyms:	1300003D20Rik; mFLJ00153
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023622
ORF Size:	2127 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218973).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023622.4 , NP_076111.2
RefSeq Size:	2550 bp
RefSeq ORF:	2130 bp
Locus ID:	71746
UniProt ID:	Q3UYI5
Cytogenetics:	9 A3
Gene Summary:	Guanine nucleotide exchange factor (GEF) for Ral-A. Potential effector of GTPase HRas and Ras-related protein M-Ras. Negatively regulates Elk-1-dependent gene induction downstream of HRas and MEK1.[UniProtKB/Swiss-Prot Function]



[View online »](#)