

Product datasheet for **MR218975L4V**

Zfp276 (NM_020497) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Zfp276
Synonyms:	AW048709; D8Ertd370e; D8Ertd377e; Znf276
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020497
ORF Size:	1842 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR218975).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_020497.2 , NP_065243.2
RefSeq Size:	4877 bp
RefSeq ORF:	1845 bp
Locus ID:	57247
UniProt ID:	Q8CE64
Cytogenetics:	8 72.09 cM



Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]