

Product datasheet for **MR218988L3V**

Ccdc136 (NM_145574) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Ccdc136
Synonyms:	492151IK06Rik; BC006583
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PSI00092)
Tag:	Myc-DDK
ACCN:	NM_145574
ORF Size:	3006 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR218988).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_145574.3 , NP_663549.3
RefSeq Size:	3836 bp
RefSeq ORF:	3009 bp
Locus ID:	232664
UniProt ID:	Q3TVA9
Cytogenetics:	6 A3.3



Gene Summary: May play a role in acrosome formation in spermatogenesis and in fertilization.
[UniProtKB/Swiss-Prot Function]