

Product datasheet for **MR218990**

Zfp385b (NM_178723) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Zfp385b
Synonyms:	B830010L13Rik; C130013B13Rik; Zfp533; Znf385b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR218990 representing NM_178723
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAACTACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAACATGGCAACCTTCTCCGGGGCTTTGAAGAAAAGGGTTAAGAACGACAGGCCTGGGGACCAGT
 TCAGCAAGGAGAAGAAGAAGATTCTCTTCTCTTTGTGAAGTGTGCAACATCCAGCTAACTCCGCTGC
 TCAGGCCCAGGTCCATTATGACGGCAATCCCATCGAAAGCGGGTCAAGCAGCTGAGCGATGGCCAGCCC
 CCACCCCCGGTCCAGGGCTCAGTGCCACTGCTGGCTGGCCCTTGCCCTTGCCCTGGGCTGGCCGAACA
 CCAGCACAGGCAGCGCATGCCACACCACTACCTTCTGCGCTCGTGGCAGACCCACCTGATGATGCA
 GCCTTCGCTGGACATCAAACCGTTTCATGTCTTTCCAGTGGACAGTAGCTCTGCCGTGGGGCTCTTTCCA
 AACTTTAACACAATGGACCCTGTGCAGAAAGCAGTCATTAACCACACATTTGGAGTCTCCATTCCCCGA
 AGAAGAAACAAGTTATTTCTGTAACTGTGTGAGCTTCGCTTTAACTCGGATAGCCAGGCCGAGGCCCA
 CTCAAAGGAAGTAAACATGCCAAGAAGGTCAAAGCTCTAGAGGCAACGAAAAATAAACCCTAAATGGTT
 CCTTCCAAGGACAGCGCAAGGCTAATCCAGCTGCTCCATCAGGCCAGGCACAGGCGACAGCTCTGACA
 AATCAGAAGACAAAGGGAAGATAAAAGCCACCAGTTCCAGTCAGCCCTCAGGCTCCGAAGGAGGCTCCTT
 TCTCTCAAATCTGGCACAACCCCTGCCACTTGGAGCCATCGCTTCTCCCTCAAGAGCACAACCGGC
 GCTCCAGGGTCTGTTGCTGAGTCAGAAGAAGAAAAGCCAAAAATTAATCTACTGCTCACTATGCAAGG
 TGGCTGTGAATCCCTGTGACAGCTGGAGGCTCACAACACAGGATCTAAACACAAGACCATGGTTGAGGC
 TCGTAATGGAGCTGGTCCAATTAAGTCCTATCCAAGACCTGGGTCAAGATTAAAGGTGCAGAACGCGAGT
 AAGGGCTCAGGACTACAGAACAAAATGTTTCATTGTGAGATCTGTGATGTTACGTGAACCTCAGAGATTC
 AACTCAAACAGCACATTTCAAGCCGAAGGCATAAAGACCGGTGGCAGGAAAACCGCTGAAGCCGAAATA
 CAGCCCTACAACAAGCTGCAGCGGAGCCCAAGCATCTTAGCTGCAAAGCTTGCAATCCAGAAAGATCTG
 ATGAAGCCTCTGGCCCTACGTTCTGTCTTACCCCTGGCGGCCGAGCCGTGCTCTGCGCTGTAC
 TGCCACCTCGGCCATCCGCTCTCTGTTCCAGGCTGCAGCCATCCCTCCTGCGCTGTGAGGCTTGCCA
 CGGGCCATCCGGGCCAGCCCGCTCCATCCTCTTTGCCCTTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218990 representing NM_178723
 Red=Cloning site Green=Tags(s)


MNMTFLRGFEKGLKNDPRGDQFSKEKKILFSFCEVCNQLNSAAQVHYDGKSHRKRKQLSDGQP
 PPPVQGSVPLLAGPCPCPGPNTSTGSACHTTLPALVRTPTLMMQPSLDIKPFMSFPVDSSSAVGLFP
 NFNTMDPVQKAVINHTFGVSIIPKKKQVISCNVCQLRFNSDSQAEAHYKGSKHAKVKALEATKNPKMV
 PSKDSAKANPSCSIRPGTGDSSDKSEDKGIKATSSSQPSGSEGSFLLKSGTTPLPLGAIASPSKSTNG
 APGSVAESEEKAKLLYCSLCKVAVNSLSQLEAHNTGSKHKTMEARNAGAPIKSYPRPGSRLKVQNGS
 KGSLQNKMFHCEICDVHVNSEIQLKQHISRRHKDRVAGKPLPKYSPYNKLQSPSILAALAFQKDL
 MKPLAPTFLSSPLAAAVSSALSLPPRPSASLFQAAIIPALLRPGHPIRATPASILFAPY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9082_d12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_178723

ORF Size: 1446 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_178723.6](#), [NP_848838.2](#)

RefSeq Size: 3217 bp

RefSeq ORF:	1449 bp
Locus ID:	241494
UniProt ID:	Q8BXJ8
Cytogenetics:	2 C3
MW:	51.4 kDa
Gene Summary:	May play a role in p53/TP53-mediated apoptosis.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MR218990