

Product datasheet for MR219003

Spata3 (NM_028647) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spata3 (NM_028647) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Spata3
Synonyms: 1700011N12Rik; 1700029H01Rik; 4930424D10Rik; Mtsarg1; Tsarg1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR219003 representing NM_028647
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAAGGTCAAAAAGAAAAGTCTGATTCCAGGCGCCGAGGAATTCATTTACCGCAGACTAGCT
CCGATTCTCACAGCAGCCAGCTCAGAAACCCCGCCGTCATGTCCAGAGCCCGCTCGCCACCATCTAA
GCCACAGCCTTGCCAGGAGTCCACCACACCACCAAGTTAACCTGAGCCTAACCGCAACAGCATACA
CCACAGCCTCTCCGCTCCAGAGAACTGCCTCATCTCGTTCCTGGTACCCATGGAACCTAAACCAA
TTCTTCAAGCAGGAAAGCAGTGCCTTGACCTATGTGGCCCCAGCTCCTGCTCCTGTGCAGCCTGCC
GGCAGCTCCGCTTGCTGGCATCGTCTGGGGCTTTGCCATAGCCGAATCTTTGATGTCCTCTACCTCGG
GACTGGTCATCCATGCCAGGGAGAGGAGTCCAAACCTTCTCACCTTCTACAGAAAACCGTCAAGGAAGT
ACTGCGCTCCCGTAACTCTCGTCTTCAAGCTCTCGAACTGTTGCTGTGGCTCTGGGGCCTTGGGAG
TTGCCTACTTCATGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219003 representing NM_028647
 Red=Cloning site Green=Tags(s)

MKKVKKKSDSRRRRNSISPQTSSDSSQQPSSETPPSCPEPASPPSKPQPCQESTTPHQVNPEPKPQQHT
PQPLPPEKPASSPFLVPMPEPKILPSRKAVPLTYVAPRSCSACPGSSACWHRLGLCHSRIFDVLLPR
DWSSMPGRGVPNLLTFYRKPSRKYCAPRNSRASSRNCCC GSGGLGSCLLHG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_028647

ORF Size: 576 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028647.5](#)

RefSeq Size: 1051 bp


RefSeq ORF: 579 bp

Locus ID: 70060

UniProt ID: [Q9D9T6](#)

Cytogenetics: 1 C5
 MW: 21.3 kDa

Product images:



Circular map for MR219003