

## Product datasheet for **MR219014L3V**

### Catsper4 (NM\_001130030) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Catsper4 (NM_001130030) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Catsper4
Synonyms:	4933412L05
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130030
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219014).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001130030.1</a> , <a href="#">NP_001123502.1</a>
RefSeq Size:	1511 bp
RefSeq ORF:	1110 bp
Locus ID:	329954
UniProt ID:	<a href="#">Q8BVN3</a>
Cytogenetics:	4 D3



[View online »](#)

**Gene Summary:**

Voltage-gated calcium channel that plays a central role in sperm cell hyperactivation. Controls calcium entry to mediate the hyperactivated motility, a step needed for sperm motility which is essential late in the preparation of sperm for fertilization. Activated by intracellular alkalinization.[UniProtKB/Swiss-Prot Function]