

Product datasheet for **MR219017L3V**

Acap2 (NM_030138) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acap2 (NM_030138) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Acap2
Synonyms:	4832442G16; 9530039J15Rik; Centb2; CNT-B2; mKIAA0041
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030138
ORF Size:	2310 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219017).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030138.2 , NP_084414.1
RefSeq Size:	6493 bp
RefSeq ORF:	2313 bp
Locus ID:	78618
UniProt ID:	Q6ZQK5
Cytogenetics:	16 B2
Gene Summary:	GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).[UniProtKB/Swiss-Prot Function]



[View online »](#)