

Product datasheet for MR219025L3

Ece2 (Eeflakmt4-ece2) (NM_177941) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Eeflakmt4-ece2
Synonyms:	1810009K13Rik; 6330509A19Rik; 9630025D12Rik; BBI27715; Ece-2; Eeflakmt4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219025).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.

ACCN:	NM_177941
ORF Size:	2730 bp



OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_177941.1 , NP_808810.1
RefSeq Size:	3436 bp
RefSeq ORF:	2733 bp
Locus ID:	110599584
UniProt ID:	PODPE0
Cytogenetics:	16 B1
Gene Summary:	This locus represents naturally occurring readthrough transcription between the adjacent genes eukaryotic translation elongation factor 1 alpha lysine specific methyltransferase 4 (GeneID:110599566) and endothelin converting enzyme 2 (GeneID:107522). The readthrough transcript representing this gene encodes a fusion protein that shares sequence identity with each individual gene product. [provided by RefSeq, Jul 2017]

Product images:



Circular map for MR219025L3