

Product datasheet for **MR219032L3V**

Myo1g (NM_178440) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Myo1g
Synonyms:	E430002D17Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178440
ORF Size:	3072 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR219032).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_178440.4 , NP_848534.2
RefSeq Size:	3307 bp
RefSeq ORF:	3075 bp
Locus ID:	246177
UniProt ID:	Q5SUA5
Cytogenetics:	11 A1



Gene Summary:

Unconventional myosin required during immune response for detection of rare antigen-presenting cells by regulating T-cell migration (PubMed:25083865). Unconventional myosins are actin-based motor molecules with ATPase activity and serve in intracellular movements. Acts as a regulator of T-cell migration by generating membrane tension, enforcing cell-intrinsic meandering search, thereby enhancing detection of rare antigens during lymph-node surveillance, enabling pathogen eradication (PubMed:25083865). Also required in B-cells, where it regulates different membrane/cytoskeleton-dependent processes (PubMed:24310084). Involved in Fc-gamma receptor (Fc-gamma-R) phagocytosis (PubMed:23038771).[UniProtKB/Swiss-Prot Function]