

Product datasheet for **MR219033L3V**

Lrrfip2 (NM_027742) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrrfip2 (NM_027742) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrrfip2
Synonyms:	5133400F20Rik; AI850587
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027742
ORF Size:	1200 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219033).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027742.3 , NP_082018.1
RefSeq Size:	3215 bp
RefSeq ORF:	1203 bp
Locus ID:	71268
UniProt ID:	Q91WK0
Cytogenetics:	9 F3



[View online »](#)

Gene Summary:

May function as activator of the canonical Wnt signaling pathway, in association with DVL3, upstream of CTNNB1/beta-catenin. Positively regulates Toll-like receptor (TLR) signaling in response to agonist probably by competing with the negative FLII regulator for MYD88-binding (By similarity).[UniProtKB/Swiss-Prot Function]