

Product datasheet for **MR219058L3V**

Ssu72 (NM_026899) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ssu72 (NM_026899) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ssu72
Synonyms:	1190002E22Rik; 1500011L16Rik; 2610101M12Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026899
ORF Size:	582 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219058).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026899.3 , NP_081175.2
RefSeq Size:	1478 bp
RefSeq ORF:	585 bp
Locus ID:	68991
UniProt ID:	Q9CY97
Cytogenetics:	4 E2
Gene Summary:	Protein phosphatase that catalyzes the dephosphorylation of the C-terminal domain of RNA polymerase II. Plays a role in RNA processing and termination. Plays a role in pre-mRNA polyadenylation via its interaction with SYMPK.[UniProtKB/Swiss-Prot Function]


[View online »](#)