

Product datasheet for **MR219078**

Armc8 (NM_028768) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Armc8 (NM_028768) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Armc8
Synonyms:	1200015K23Rik; Gid5; HSPC056
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR219078 representing NM_028768
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCGTGCCTGTTGGAGACCCCAATCCGCATGAGCGTCTCTCGGAGGTAACAGCTAGCAGTCCCACT
ATGTTGACAGGCTATTTGACCCTGATCCCCAGAAAGTTCTCCAAGGTGTCATAGACATGAAAAATGCTGT
AATTGGAAACAATAAACAGAAAGCCAACCTCATCGTTCTAGGAGCGGTTCCAAGATTGTTGTATTTGCTT
CAACAAGAAACCTCAAGTACAGAGCTGAAAAGTGAATGTGCAGTGGTGTAGGGAGTCTTGCGATGGGCA
CAGAAAACAATGTCAAGTCTCTACTAGACTGCCATATTATCCCTGCCTTACTACAAGGACTACTGTCCCC
AGACCTGAAGTTCATTGAAGCTTGCCCTCCGATGCCTCCGACAATCTTACCAGCCCTGCTACTCCAGAG
GAGCTACTGTATACAGATGCTACGGTGATACCACACCTCATGGCACTGCTGAGCAGGTCCCGCTACACCC
AGGAGTATATCTGCCAGATCTTCTCACACTGTTGCAAAGGGCCAGATCATCAAACAATTTTGTAAACCA
TGGGGCAGTGCAGAACATTGCTCATCTCCTAACCTCACCTCTTACAAAGTTCGAATGCAAGCACTGAAA
TGCTTCTCTGTTTTAGCCTTTGAGAATCCCAAGTATCGATGACCCTGGTAAACGTTTTGGTTGATGGGG
AGTTGTTACCACAGATATTTGTGAAGATGTTACAGAGGGATAAGCCTATTGAGATGCAGCTCACGTCAGC
AAAATGTTAACTTACATGTGCAGAGCTGGAGCAATTCGCACAGATGACAGCTGTATTGTATAAAACT
TTGCCCTGTTTGGTTTCAATGTGCAGTAAGGAGAGGTTACTAGAGGAGAGAGTTGAAGGAGCTGAGACCC
TCGCCTATCTGATTGAGCCAGATGTTGAGCTGCAGAGAATCGCTAGCATAACCGACCACCTCATTGCCAT
GCTGGCTGACTATTTCAAGTATCCCAGCTCAGTGCAGGCCATCACTGACATTAAGGCTTGATCATGAC
TAAAAACATGCCCATGAGCTCCGCCAGGCTGCGTTCAAGCTCTATGCCTCTCTGGAGCAAAATGATGAAG
ACATCCGGAAGAAGATCATTGAGACTGAAACTATGATGGACCGAATTGTGACTGGCTTGTCTGAGTCTAG
TGTCAAGGTGCGGTTAGCTGCTGTCAGATGTTTGCACAGTTTATCCAGATCTGTACAACAGCTTCGAACC
AGTTTCCAGGATCATGCGGTGTGAAAACCTTTAATGAAGTTTTACAAAATGCACCCGATGAAATCCTAG
TAGTAGCATCTTCCATGTTATGTAACCTTCTTCTGAATTTTCTCCAAGTAAAGACCAATTTTGGAAATC
AGGAGCTGTCGAGCTGCTCTGTGGGTTGACACAGAGTAAAAACCTGCTCTACGAGTGAATGGGATTTGG
GCTTTAATGAATATGGCATTTCAGGCGGAACAAAAATAAAGCAGATATTTACGAAGTCTGAGTACTG
AACAGTTATTCGGTTATTGTCAGATTCTGATATGAATGTGCTGATGAAGACATTGGGGCTTCTTAGAAA
TCTCCTCTCCACTCGGCCTCATATAGATAAAATCATGAGTACTCATGGGAAGCAGATTATGCAAGCTGTC
ACTCTGATACTGGAAGGGGAGCATAGCATTGAAGTCAAAGAGCAGACCCTGTGCATCCTAGCCAACATTG
CAGATGGGACAACAGCAAAGAGCTTATTATGACCAATGATGACATCCTGCAGAAAATCAAATATTACAT
GGGTCAATTCACATGTCAAACGAGCTTGTGCCATGTTTTGTATATCAAACCTCATATGGAACGAAGAG
GAAGGTTACAGGAACGCCAGGATAAATTACGAGACATGGGGATTGTAGATATTCTACACAAACTGAGTC
AGTCGGCAGACTCAAACCTCTGTGACAAGGCAAAGACGGCGCTGCAGCAGTACCTGGCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219078 representing NM_028768
Red=Cloning site Green=Tags(s)

```

MACLLETPIRMSVLSEVTASSRHYVDRLFDPDPQKVLQGVIDMKNNAVIGNNKQKANLIVLGA VPRLLYLL
QQETSSSTELKTECAVVLGSLAMGTENNVKSLLDCHIIPALLQGLLSPDLKFEIACLRLCRTIFTSPVTP
ELLYTDATVIPHLMALLSRRYTQEYICQIFSHCCKGPDHQTILFNHGAVQNIHLLTSPSYKVRMQALK
CFSVLAFENPQVSM TLVNVLV D GELL PQIFVKMLQRDKPIEMQLTS AKCLTYMCRAGAIRTD DDCIVLKT
LPCLVRMCSKERLLEERVEGAETLAYLIEPDVELQRIASITDHLIAMLADYFKYPSSVSAITDIKRLDHD
LKHAHEL RQA AFKLYASLGANDEDIRKKI IETETMMDRIVTGLSESSVKVRLAAVRCLHSLSRVQQRLRT
SFQDHAVWKPLMKVLQNAPDEILVVASSMLCNLLLEFSPSKEPILESGAVELLCGLTQSEN PALRVNGIW
ALMNMFAQAEQKIKADILRSLSTEQLFRLLSDSDMNVLMKTLGLLRNLLSTRPHIDKIMSTHGKIQM QAV
TLILEGEHSIEVKEQTLCILANIADGTTAKELIMTND DILQKIKIYYMGHSHVKLQLAAMFCISNL IWNEE
EGSQERQDKLRDMGIVDILHKLQSADS NLCDKAKTALQQYLA
    
```

TRTRPLEQKLISEEDLANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9090_h03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_028768

ORF Size: 2019 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028768.3](#), [NP_083044.2](#)

RefSeq Size: 4628 bp

RefSeq ORF: 2022 bp

Locus ID: 74125


UniProt ID: [Q9DBR3](#)

Cytogenetics: 9 E3.3

MW: 75.4 kDa

Gene Summary: Component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219078