

## Product datasheet for **MR219078L4V**

### Armc8 (NM\_028768) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Armc8 (NM_028768) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Armc8
Synonyms:	1200015K23Rik; Gid5; HSPC056
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028768
ORF Size:	2019 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219078).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028768.3</a> , <a href="#">NP_083044.2</a>
RefSeq Size:	4628 bp
RefSeq ORF:	2022 bp
Locus ID:	74125
UniProt ID:	<a href="#">Q9DBR3</a>
Cytogenetics:	9 E3.3
Gene Summary:	Component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1.[UniProtKB/Swiss-Prot Function]



[View online »](#)