

## Product datasheet for **MR219082L4V**

### **Snrpb (NM\_009225) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Snrpb (NM_009225) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snrpb
Synonyms:	AL024368; AU018828; SM-B; SM11; SMB; SNRNP-B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009225
ORF Size:	693 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219082).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009225.2</a> , <a href="#">NP_033251.1</a>
RefSeq Size:	1121 bp
RefSeq ORF:	696 bp
Locus ID:	20638
UniProt ID:	<a href="#">P27048</a>
Cytogenetics:	2 63.19 cM



[View online »](#)

**Gene Summary:**

Plays role in pre-mRNA splicing as core component of the SMN-Sm complex that mediates spliceosomal snRNP assembly and as component of the spliceosomal U1, U2, U4 and U5 small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Component of both the pre-catalytic spliceosome B complex and activated spliceosome C complexes. Is also a component of the minor U12 spliceosome (By similarity). As part of the U7 snRNP it is involved in histone pre-mRNA 3'-end processing (PubMed:19470752). [UniProtKB/Swiss-Prot Function]