

Product datasheet for **MR219085**

Krt84 (NM_008474) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Krt84 (NM_008474) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Krt84
Synonyms:	AV235125; HRb-1; Krt2-3; Krt2-16
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR219085 representing NM_008474
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCCTTGCCGCTCCTACAGAGTCAGCTCTGGCCGTCGGGTGGGAAGCTTCAGTTCCTGCTCAGCAATGG
 CCCCTCAGCATCTGAACCGCTTCCGGTCCAGCTCTGTCTCCTGCAGGAGTGGCCTGGCTTCCGGGGCTT
 GGGTGGTTTTGGTAGTCGGAGCGTCATCAACTTTGGATCATCCTACCCCGGATAGCAGTTGGATGCTCT
 CGTCCTGTCCGCTGTGGAGTTGGCTTTGGAGCTGGCAGTGGGATGGCCTTTGGATCTGGTATGGCCTGG
 GCTTTAGGGCCAGCAGCGGTGTTGGCCTGGGTTTGGGCTGGCGGCTGCCCTAGCTATGGATTCGGAGG
 CCCTGGATTGGAGGCCCTGGCTTTGGTGGCCCTGGCTTTGGCTACAGAATCGGAGGAATTGGAGGCCCA
 TCAGCCCCGCTATCACAGCTGTAAGTTAAACCAGAGCCTGCTGACACCCCTCAACCTGGAGATTGACC
 CCAATGCCAGAGGGTGAAGAGAGACGAAAAGGAGCAGATCAAGACCCTCAACAACAAATTTGCCTCCTT
 CATTGATAAGGTGAGTTTCTGGAGCAACAGAATAAACTCCTAGAGACCAAGTGGAGCTTCTCCAAGAA
 CAGAAATGTGCCCGGAGCAATCTGGAGCCTCTCTTCGATAACTATATACCAACCTGCGGAGGCAGCTGG
 ATGTCCTAAGCAGTGACCAAGCTCGGCTGCAGGCAGAGAGGAACCACTGCAAGACATCCTTGAGGGCTT
 CAAGAAGAAGTATGAAGAGGAGGTGATTCCGGGCCAATGCCGAGAATGAGTTTGTAGCGCTGAAGAAG
 GATGTAGATGCAGCTTCTCTGAACAAATCTGATCTGGAAGCCAATGTGGATACTCTGATCCAGGAACTG
 AATTCCTGAAGGCTCTCTACCATGAGGAAATCGAGATGCTGCAATCACACATATCAGAGACATCTGTCAT
 CGTGAAGATGGACAACAGCCGGGATCTGAACCTAGATGGGATCATAGCTGAAGTCAAGGCACAGTATGAG
 GAGGTGGCGAGGCGGAGTCGTGCAGATGTGGAGTCGTGGTACCAGACCAAGTATGAGGAGATCGGGTGA
 CCGCTGGCCAGCACTGTGACAATCTACGTAGCACACGGGACGAGATCAATGAAGTACCCGGTTGATCCA
 GAGGCTGAAAGCAGAGATCGAGCATACCAAGGCTCAGTGTGCCAACTGGAGGCTGCTGTGGCTGAGGCA
 GAGCAGCAGGGTGAAGCAGCCCTCAATGATGCCAAGTGAAGCTGGCTGATCTGGAGGGTGCCTGCAGC
 AGGCCAAGCAGGACATGGCCAGGCAGCTGCGGGAGTACCAGGAGCTCATGAATGTCAAGCTGGCCCTGGA
 CATCGAGATTGTACCTACAGGAGGCTGCTGGAGGCGAGGAAATTCGGATCTGCGAAGGAGTTGGGCCA
 GTAAACATATCTGTGAGCAGCTCTCGAGGAGGTGTGCTGTGTGGGCCGAGTTAGTCTCTGGATCCAGCC
 TTTCCACAATGGTGGAGTACCTTCTCCACCAGCAGCAGCATCCGTGCCACCGGAGGGGTTCTGGCTTC
 CTCAAGCCTGAGGGCTGGTGGGACCTGCTGAGCTCTGGGTCCAGAGGAGGCTCAGTGTGGTGGGTGAT
 GCCTGTGCCCCCAGCATCCCCTGTGCACTGCCCCACTGAGGAGGCTTCAGTAGCTGCAGTGGTGGTCTG
 GTAACAGGAGCTCCAGCTCCGCTTCTCATCTACCACCACCTCCCGCAGGACCAGGTAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR219085 representing NM_008474
 Red=Cloning site Green=Tags(s)

MSCRSYRVSSGRRVGSFSSCSAMAPQHLNFRSSSVSCRSGPGFRGLGGFGSRSVINFGSSSPRIAVGCS
 RPVRCGVGFAGSGMAFGSGDGLFRASSGVGLFGAGGCPSTYFGGPGFGGPGFGGPGFYRIGGIGGP
 SAPAITAVTVNQSLLTPLNLEIDPNAQVRKRDEKEQIKTLNNKFASFIDKVRFLQKQNKLETKWSFLQE
 QKCARSNLEPLFDNYITNLRQLDVLSSDQARLQAERNHLQDILEGFKKKYEEVVFRANAENEFVALKK
 DVDA AFLNKS DLEANVDTLIQETEFKALYHEEIEMLQSHISETSVIVKMDNSRDLNLDGIIAEVKAQYE
 EVARRSRADVSWYQTKYEEMRVTAGQHCDNLRSTRDEINELTRLIQRLKAEIEHTKAQCAKLEAAVAEA
 EQQGEAALNDAKCKLADLEGALQQAQDMARQLREYQELMNVKLALDIEIVTYRRLLEGEIEIRICEGVGP
 VNISVSSSRGGVLCGPVLVSGSSLSHNGGVTFSTSSSIRATGGVLASSSLRAGGDLLSSGSRGGSVLVGD
 ACAPSIPCALPTEGGFSSCSGGRGNRSSSVRFSSSTTSRRTRY


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mm9044_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_008474

ORF Size: 1809 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_008474.2](#), [NP_032500.2](#)

RefSeq Size: 2570 bp

RefSeq ORF: 1812 bp

Locus ID: 16680
 UniProt ID: [Q99M73](#)
 Cytogenetics: 15 56.9 cM
 MW: 65.4 kDa

Product images:



Circular map for MR219085