

## Product datasheet for MR219127

### Lamp5 (NM\_029530) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Lamp5 (NM_029530) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lamp5
Synonyms:	3110035N03Rik; 6330527O06Rik; BAD-LAMP; LAMP-5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR219127 representing NM_029530 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGATCTCCGAGTAAGAACCCTTCTTGGTGGAGACAGACTTCGGATTCTCCTAATGTTCTTCCATGTAA  
TGGTCCAAACCGTGCCAGAACAAGAAGTGGAAAATCTTTCGGGCCTTCCACGAACCCCGAGAAAGACAT  
ATTTGGTGGTGGGAAAATGGGACGACGTGTCTCATGGCAGAGTTGCAGCCAAATTTATTGTACCTTAT  
GATGTGTGGGCCAGCAATTATGTGGATCTGATCACAGAACAGGCTGAGATCTCCTTGACCCGGGGAGCTG  
AGGTGAAGGGCCACTGTGGCCACAACGAGTCGGAGCTGGAGGTGTCTGGGTAGATCACGCCTACACACT  
CAGAATGCTCTTTGTAAGGAAAAGTCACAACACTTCCAAAGGACCAGAGGCGACTTGGAACTGAACAAG  
GTGCATTTTGTCTATGATTCCTCGGAGAAGACTCACTTTAAAGCCCTGTGAAGGTTAATAAGTACATCG  
CCAGCTCTCATCACCTCTCTGCCTTGGTCAACCCAGCTGGGATGTCTATGAGTGTGAGGCTCAGCAGAC  
CATTTCTCTGGCCTCCAGTGACCCTCAGAAGACTGTCACCATGATCCTGTCCGAGTGCACATCCAGCCC  
TTTGATATCATCTCTGACTTTGTCTTCAAGTGAAGAACAATAATGTCCAGTGGATGAGCAAGAGCAGTTGG  
AAGAGACCCTGCCCTGATCCTGGGCCTCATCTGGGCCTTGTATTGTCATAACACTTGTGATCTACCA  
CATCCACCATAAAATGACTGCCAATCAAGTGCAGATCCCAGGGACCGTTCAGTACAAGCACATGGGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR219127 representing NM\_029530  
 Red=Cloning site Green=Tags(s)

MDLRVRTLLGGDRLRILLMFFHVMVQTVAEQEVENLSGLSTNPEKDI FVVRENGTTCMAEFAAKFIVPY  
 DVWASNYVDLITEQAEISLTRGAEVKKGHCGHNESELEFVVDHAYTLRMLFVKESHNTSKGPEATWNLNK  
 VHFVYDSSEKTHFKAPVKVNKYIASSHLSALVTPAGMSYECQAQQTISLASSDPQKTVTMILSAVHIQP  
 FDIISDFVFSEEHKCPVDEQEQL EETLPLILGLILGLVIVITLVIYHIHHKMTANQVQIPRDRSQYKHMG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9029\\_a11.zip](https://cdn.origene.com/chromatograms/mm9029_a11.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_029530

**ORF Size:** 840 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_029530.3](#)

**RefSeq Size:** 1965 bp

**RefSeq ORF:** 843 bp

**Locus ID:** 76161


**UniProt ID:** [Q9D387](#)

**Cytogenetics:** 2 F3

**MW:** 32.2 kDa

**Gene Summary:** Plays a role in short-term synaptic plasticity in a subset of GABAergic neurons in the brain. [UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR219127