

## Product datasheet for **MR219138L4V**

### Ankrd16 (NM\_177268) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ankrd16 (NM_177268) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ankrd16
Synonyms:	2810455F06Rik; A1646698; D430029B21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177268
ORF Size:	1083 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219138).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_177268.4</a> , <a href="#">NP_796242.2</a>
RefSeq Size:	2356 bp
RefSeq ORF:	1086 bp
Locus ID:	320816
UniProt ID:	<a href="#">A2AS55</a>
Cytogenetics:	2 A1



[View online »](#)

**Gene Summary:**

Required to prevent the misactivation of serine (Ser) with tRNA(Ala) by promoting the hydrolysis of Ser-mischarged tRNA(Ala), thereby playing a role in translational fidelity (PubMed:29769718). Binds directly to the catalytic domain of AARS/AlaRS and captures Ser that is misactivated by AARS/AlaRS, preventing the charging of Ser adenylates to tRNA(Ala) and precluding Ser misincorporation in nascent peptides (PubMed:29769718). [UniProtKB/Swiss-Prot Function]