

Product datasheet for **MR219143L3V**

Mex3c (NM_001039214) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mex3c (NM_001039214) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mex3c
Synonyms:	A130001D14Rik; BC035207; BM-013; Rkhd2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039214
ORF Size:	1956 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219143).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039214.4 , NP_001034303.3
RefSeq Size:	4080 bp
RefSeq ORF:	1959 bp
Locus ID:	240396
UniProt ID:	Q05A36
Cytogenetics:	18 E2



[View online »](#)

Gene Summary:

RNA-binding protein. May be involved in post-transcriptional regulatory mechanisms, modulating levels of some mRNAs by promoting their degradation in a way involving ubiquitin ligase activity. May act as suppressor of replication stress and chromosome missegregation.[UniProtKB/Swiss-Prot Function]