

Product datasheet for **MR219145L3V**

Ubap2l (NM_001165986) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ubap2l (NM_001165986) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ubap2l
Synonyms:	3110083O19Rik; 4932431F02Rik; A430103N23Rik; C77168; mKIAA0144; Nice-4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001165986
ORF Size:	3042 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219145).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001165986.1 , NP_001159458.1
RefSeq Size:	3720 bp
RefSeq ORF:	3045 bp
Locus ID:	74383
UniProt ID:	Q80X50
Cytogenetics:	3
Gene Summary:	Plays an important role in the activity of long-term repopulating hematopoietic stem cells (LT-HSCs) (PubMed:25185265).[UniProtKB/Swiss-Prot Function]



[View online »](#)