

Product datasheet for **MR219150L4V**

Saa2 (NM_011314) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Saa2 (NM_011314) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Saa2
Synonyms:	AW111173; Saa-2; Saa1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011314
ORF Size:	366 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219150).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011314.2
RefSeq Size:	610 bp
RefSeq ORF:	369 bp
Locus ID:	20209
UniProt ID:	P05367
Cytogenetics:	7 30.56 cM
Gene Summary:	Major acute phase reactant. Apolipoprotein of the HDL complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)