

Product datasheet for **MR219156L3V**

Acrbp (NM_016845) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acrbp (NM_016845) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Acrbp
Synonyms:	OY-TES-1; sp32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016845
ORF Size:	1620 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219156).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016845.2 , NP_058541.2
RefSeq Size:	1913 bp
RefSeq ORF:	1623 bp
Locus ID:	54137
UniProt ID:	Q3V140
Cytogenetics:	6 F2
Gene Summary:	May be involved in packaging and condensation of the acrosin zymogen in the acrosomal matrix via its association with proacrosin.[UniProtKB/Swiss-Prot Function]



[View online »](#)