

Product datasheet for **MR219174**

Gpr160 (NM_027965) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gpr160 (NM_027965) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gpr160
Synonyms: 1700025D19Rik; GPCR150
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR219174 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAGCTCTTTCTTCAAAGAAGTCTCTCTTCAATACCAGCTCCATCAGTCTCCGAGCTCCTGGAGG
CCTCCTGTCTGTTCTTGATCATCTGGGAAAGTGTATTAAACATCCTCCTCTAAGAGTGAGAAG
GGGAGACGCCCGCTGGACTTTGATGGAGTACTTCTGCTTCTCACTCGCACTTGTCGATCTCCTGCTGTTG
GTAACATTTCCATCCTGACCTACTTCAGGGACTTTGTAGTTCTAGGTATTAGGTTTACCAGATACCACA
TCTGTCTGCTCACACAAATTATTTCTTCACTTATGGCTTCTTGCAATTACCCAGTCTGCTCACTGGCTTG
TATAGACTACTGGTGAATCTGTCCAGAGCCAGCAAGCAGTCATCTAGGTGGCAGAAATTAATTTATTTT
TTGACAGTCATTTAACTTGATCTCAGTCCTTGCTTACGTTTTGGTTGACCCGGCCATCTCAGTCAGCC
TAAAGGCACACAGGGGCTATGTGTACCAGTGCCCGCCTATGTCAGCACTCAGAGTCACTGGCTCTCACT
GTCTATGCTGATGGTTTTATTTGTGGCTTTTCTGATTTTCATGGCAGGAAGTTGTTGCCTTGCTTCAGGCT
ATGAGGATAGCATCCTATAAGAGCAAGGCCGCTCTACTTTCTTTCCCACTCCACTGCGGTTACGCC
TGAGCTGCAGAGAAGCACTCTTGCCAGGCTCATCGTGTGCTTTCTTGGACCTGGTTTCCCTTTGTAGC
ACTTCAGTCCATCCTCTCCCTCAGAGTTTTCAGATCCCGGCGTACATAGAAATGAACGTCCCCTGGCTG
TACTTTGTCAACAGTTTCTCATCGCTGCAGTTTATTGGTTTAACTGCCACAAGTTGATTTAAGAGACA
GCTCGTTACCTGTGGATCCTTTCATCAACTGAAATGCTGCTTTGTTCCAGTCCATAGACTTAAGCAAGT
GGAGAGGCCCATGTCCATAGTCATCTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR219174 protein sequence
Red=Cloning site Green=Tags(s)

MTALSSKNCSLQYLHQSPQLLEASCLLFLIILGKVLNILLRVRRGDARWTLMEYFCFSLALVDLLLL
 VNISILTYFRDFVVLGIRFTRYHICLLTQIISFTYGFLHYPVCSLACIDYWCNLSRASKQSSRWQKLLYF
 LTVILTWISVLAYVLVDPAISVSLKAHRGYVYQCPAYVSTQSHWLSLSMLMVLVFAVLISWQEVVALLQA
 MRIASYKSKAALYFPFPLHCGYALSCEALLPRLIVCFLGTWFFPVALQVLILSLRVQIPAYIEMNVPWL
 YFVNSFLIAAVYWFNCHKLDLDRDSSLPVDPFINWKCCFVPVHRLKQVERPMSIVIC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_027965

ORF Size: 1008 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_027965.2](#), [NP_082241.1](#)

RefSeq Size: 2020 bp

RefSeq ORF: 1011 bp

Locus ID: 71862


UniProt ID: [Q3U3F9](#)

Cytogenetics: 3 A3

MW: 38.9 kDa

Gene Summary: Orphan receptor.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219174