

## Product datasheet for **MR219177**

### **Tax1bp3 (NM\_029564) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tax1bp3 (NM\_029564) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tax1bp3  
**Synonyms:** 1300011C24Rik; TIP-1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR219177 ORF sequence  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCTACACCCGGGCCAGCCTGTCACCGCCGTAGTGCAAAGAGTTGAAATTCATAAGTTGCGTCAAG  
GTGAGAACTTAATCTTGGGCTTCAGTATTGGAGGTGGGATCGACCAGGACCCGTCTCAGAATCCCTTCTC  
GGAAGATAAAACAGACAAGGGCATTACGTCACACGAGTATCAGAGGGAGGTCCTGCTGAAATTGCTGGG  
CTGCAGATTGGAGACAAGATCATGCAGGTGAATGGCTGGGACATGACCATGGTCACTCACGACCAGGCTC  
GGAAGCGGCTACCAAGCGCTCGGAGGAGTGGTCCGCCTGCTGGTGACTCGGCAGTCTCTACAAAAGGC  
TGTACAGCAGTCCATGCTGTCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR219177 protein sequence  
**Red=Cloning site Green=Tags(s)**

MSYTPGQPVTAVVQRVEIHKLRQGENLILGFSIGGGIDQDPSQNPFSSEKTDKGIYVTRVSEGGPAEIAIG  
LQIGDKIMQVNGWDMTMVTHDQARKRLTKRSEEVVRLLVTRQSLQKAVQQSMLS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_029564

**ORF Size:** 375 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_029564.2](#), [NP\\_083840.1](#)


**RefSeq Size:** 1417 bp

**RefSeq ORF:** 375 bp

**Locus ID:** 76281

**UniProt ID:** [Q9DBG9](#)  
**Cytogenetics:** 11 B4  
**MW:** 13.7 kDa  
**Gene Summary:** May regulate a number of protein-protein interactions by competing for PDZ domain binding sites. Binds CTNNB1 and may thereby act as an inhibitor of the Wnt signaling pathway. Competes with LIN7A for KCNJ4 binding, and thereby promotes KCNJ4 internalization. May play a role in the Rho signaling pathway (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MR219177