

Product datasheet for **MR219240L4V**

P2rx6 (NM_001159561) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2rx6 (NM_001159561) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	P2rx6
Synonyms:	P2rxl1; P2x6; P2xm
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159561
ORF Size:	1086 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219240).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001159561.1 , NP_001153033.1
RefSeq Size:	2281 bp
RefSeq ORF:	1089 bp
Locus ID:	18440
UniProt ID:	O54803
Cytogenetics:	16 A3
Gene Summary:	Receptor for ATP that acts as a ligand-gated ion channel.[UniProtKB/Swiss-Prot Function]



[View online »](#)