

## Product datasheet for **MR219317L3V**


### Slc43a3 (NM\_021398) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc43a3 (NM_021398) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc43a3
Synonyms:	Eeg1; SEEEG-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021398
ORF Size:	1509 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219317).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021398.2</a> , <a href="#">NP_067373.2</a>
RefSeq Size:	2601 bp
RefSeq ORF:	1509 bp
Locus ID:	58207
UniProt ID:	<a href="#">A2AVZ9</a>
Cytogenetics:	2 D
Gene Summary:	Putative transporter.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Product images:**

[MR219317L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR219317L3V particle to overexpress human Slc43a3-Myc-DDK fusion protein.