

Product datasheet for **MR219361**

Slc6a14 (NM_020049) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc6a14 (NM_020049) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc6a14
Synonyms:	1110007A17Rik; 9030613J17Rik; ATB0plus; CATB0plus
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR219361 representing NM_020049
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGACAGATTGAAGTGCCGAACCTTCTCAAGTGCAGACAGAAGGAGAAAGTGACAGCTTCATCCGAGA
ACTTCCATGTTGGTGAAAATGATGAAAATCAGGAACGTGGTAACTGGTCCAAAAAGTCAGATTATCTTCT
TTCATGTTGGCTATGCAGTGGGATTGGGCAATGTGTGGAGATCCCATATCTGACCTACACCAATGGC
GGAGGTGCTTTCTTGATACCCTATGCAATCATGCTAGCATTGGCTGGTCTGCCTTTGTTCTTTCTGGAAT
GTTCACTGGGACAATTTGCTAGCTTAGGTCCAGTTTCAGTTTGGAGATTTTGCCACTGTTTCAAGGTGT
GGGAATCACGATGGTCTGATATCCGTTTTGTGGCAATCTATTATAATGCATAATTGCCTACAGCCTT
TACTACTGTTTCTTTTCAAAGTGTATTGCCGTGGGCAAATGTTCTTCTGGGCAGACGAAAACCT
GTAGCAGAACACCAATAGTAACCTGGCTGTAATGTGAGCATAGGGGCAGGAGAGATGTTTATGAATATAAG
CTGGGTAAATACCAACAATCTCACCTGCTGAATGGTAGTGAAGTTTTTCGGCCAGGACAACCTCCAGT
GAACAGTATTGGGATAAAGTGACACTTCAGAGGTCAAGTGAATGGATGAAACTGGAGTTGTTGTGTGGT
ATTTAGCACTTTGCCTTCTTCTGGCTTGGCTCATAGTTGGAGCAGCACTGTTTAAAGGAATCAAATCCTC
TGGCAAGGTGGTATATTTACAGCTCTTTCCCATATGTGGTCTTCTCATTCTGTTAATACGAGGTGCA
ACTTTGGAGGGTGCATCAAAAGGCATTTTCAATACATTGGAGCAGCACTCAAAATTTACAAAACCTAAGAG
AAGCTGAGGTCTGGAAGGATGCTGCCACTCAGATATTTTATCCCTTCTGTGGCTTGGGGTGGCCTAGT
GGCTCTGCATCTTACAACAAGTTCAACAACAAGTGTACTCTGATGCCATTATAGTTTGTCTGACAAAAC
TGCTCACTAGTGTGTTGCTGGATTTGCTATATTTTCTATTCTGGGACACATGGCTCATATATCTGGAA
AAGAAGTCTCTCAAGTTGAAAATCAGGATTTGACTTGGCATTTCCTATCCTGAAGCTCTAGCTCA
ACTCCCAGCTGGTCCATTCTGGTCCATATTATTTTCTCATGCTTTTAAACATTGGGCCTTGATTCTCAG
TTTGCTTCTATTGAAACAATTACAACAACATTTCAAGATTTGTTCCCAAAGCAATGAAACGCATGAGGG
TTCCTATAACTTTGGGATGCTGCTTGTATTTGTTTCTCCTTGGCCTCCTCTGTGTGACTCAGGCTGGAAT
TACTGGGTTCAATTTGATTGATCACTTCTGTGCTGGATGGGGAATTTAATTGCTGCTATCTTGGAAAT
GCAGGAATCATCTGGATTTATGGAGGAAATAGATTCATTGAAGATATAGAAATGATGATTGGAGCAAAAA
GATGGATATTCTGGCTTTGGTGGAGAGCTTGTGGTTTGTGCTCACTCCTATTCTGTTGTCTGCAATTCT
AGTCTGGTCACTGGTGAATTTACAGACCTGATTATGCTGATATCCATATCCTGACTGGGGAGTTGCT
CTTGGCTGGTGTATGATCATTCTGCATTATTTGGATACCAATTATGGCTATCATAAAAAATAGTTCAGG
CAGAAGGAAACATACTTCAGCGTATCATAAGCTGCTGTAGACCAGCTTCGAACTGGGGTCCATACCTGGA
AAAACATCGTGGGAGAGATACAGAGACATGGCTGAACCTGCAAAGGAGACTGACCACGAAATACCTACT
ATTAGTGGCAGTACAAAACCAGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219361 representing NM_020049
 Red=Cloning site Green=Tags(s)


```
MDRLKCPNFFKCRQKEKVTASSENFHVGENDENQERGNWSKSDYLLSMVGYAVGLGNVWRFPPLYTYTNG
GGAFLLIPYAIMLALAGLPLFFLECSLGGFASLGPVSVWRILPLFQGVGITMVLISVFAIYYNVIYAISL
YYLFAFQSVLPWANCSSWADENCSTRPIVTGCNVSIGAGEMFMNISWVNTNNTCLNGSEVFRPGQLPS
EQYWDKVTLQRSSGMDETGVVVWYLALCLLLAWLIVGAALFKGIKSSGKVVYFTALFPYVLLILLIRGA
TLEGASKGISYYIGAQSNTKLREAEVWKAATQIFYSLSVAVWGLVALSSYNKFNNNCYSDAIIVCLTN
CLTSVFAGFAIFSILGHMAHISGKEVSQVVKSGFDLAFIAYPEALAQLPAGPFWWSILFFMLLTLGLDSQ
FASLETITTTFDLFPKAMKRMVPITLGCCLILFLLGLLCVTQAGIYWVHLIDHFCAGWGILIAAILEI
AGIIWIYGGNRFIEDIEMMIGAKRWIFLWWRACWFVITPILLSAILVWSLVKFRPDYADIPYPDWGVA
LWCMIIFCIIWIPIMAIKIVQAEGNILQRIISCCRPASNWGPYLEKHRGERYRDMAPAKETDHEIPT
ISGSTKPE
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9014_c07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_020049

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020049.4](#), [NP_064433.3](#)

RefSeq Size: 3433 bp

RefSeq ORF: 1917 bp

Locus ID: 56774


UniProt ID: [Q9JMA9](#)

Cytogenetics: X A2

MW: 71.9 kDa

Gene Summary: Mediates the uptake of a broad range of neutral and cationic amino acids (with the exception of proline) in a Na(+)/Cl(-)-dependent manner.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219361