

Product datasheet for **MR219361L4V**

Slc6a14 (NM_020049) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc6a14 (NM_020049) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc6a14
Synonyms:	1110007A17Rik; 9030613J17Rik; ATB0plus; CATB0plus
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020049
ORF Size:	1914 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219361).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020049.4 , NP_064433.3
RefSeq Size:	3433 bp
RefSeq ORF:	1917 bp
Locus ID:	56774
UniProt ID:	Q9JMA9
Cytogenetics:	X A2
Gene Summary:	Mediates the uptake of a broad range of neutral and cationic amino acids (with the exception of proline) in a Na(+)/Cl(-)-dependent manner.[UniProtKB/Swiss-Prot Function]



[View online »](#)