

## Product datasheet for **MR219368L3V**

### **Cnot4 (NM\_001164414) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cnot4 (NM_001164414) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cnot4
Synonyms:	Not4; Not4h; Not4hp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164414
ORF Size:	1926 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219368).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164414.1</a> , <a href="#">NP_001157886.1</a>
RefSeq Size:	3344 bp
RefSeq ORF:	1929 bp
Locus ID:	53621
Cytogenetics:	6 B1
Gene Summary:	Has E3 ubiquitin ligase activity, promoting ubiquitination and degradation of target proteins. Involved in activation of the JAK/STAT pathway. Catalyzes ubiquitination of methylated RBM15.[UniProtKB/Swiss-Prot Function]



[View online »](#)