

Product datasheet for **MR219380**

Tnfsf8 (NM_009403) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tnfsf8 (NM_009403) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tnfsf8
Synonyms:	Cd30l; CD30LG; CD153; Tnl3a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR219380 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCCAGGGCTGCAACAAGCAGGCAGCTGTGGGGCTCCTTCCCCTGACCCAGCCATGCAGGTGCAGC
CCGGCTCGGTAGCCAGCCCCTGGAGAAGCACGAGGCCCTGGAGAAGCACAAAGTCGCAGCTACTTCTACCT
CAGCACCACCGCACTGGTGTGCCTTGTGTGGCAGTGGCGATCATTCTGGTACTGGTAGTCCAGAAAAAG
GACTCCACTCAAATACAACAGAGAAGGCCCCCTTAAAGGAGGAAATTGCTCAGAGGATCTCTTCTGTA
CCCTGAAAAGTACTCCATCCAAGAAGTCATGGGCCTACCTCCAAGTGTCAAAGCATCTCAACAATACCAA
ACTGTCATGGAACGAAGATGGCACCATCCACGGACTCATATACCAGGACGGGAACCTGATAGTCCAATTC
CCTGGCTTGTACTTCATCGTTTGCCAACCTGCAGTTCCTCGTGCAGTGTCAAATCATTCTGTGGACCTGA
CATTGCAGCTCCTCATCAATCCAAGATCAAAAAGCAGACGTTGGTAACAGTGTGTGAGTCTGGAGTTCA
GAGTAAGAACATCTACCAGAATCTCTCTCAGTTTTTGTGTCATTACTTACAGGTCAAATCTACCATATCA
GTCAGGGTGGATAATTTCCAGTATGTGGATACAAACACTTCCCTCTTGATAATGTGCTATCCGCTCTTCT
TATATAGTAGCTCAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR219380 protein sequence
 Red=Cloning site Green=Tags(s)

MEPGLQQAGSCGAPSPDPAMQVQPGSVASPWSTRPWRSTSRSYFYLTSTALVCLVVAVAILVLVVQKK
 DSTPNNTTEKAPLKGNCSEDLFCTLKSTPSKKSWAYLQVSKHLNNTKLSWNEDGTIHGLIYQDGNLIVQF
 PGLYFIVCQLQLVQCSNHSVDLTLQLLINSKIKKQTLVTVCESGVQSKNIYQNLSQLLLHYLQVNSTIS
 VRVDNFQYVDTNTFFPLDNVLSVFLYSSSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_009403

ORF Size: 720 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009403.3](#), [NP_033429.1](#)

RefSeq Size: 2228 bp

RefSeq ORF: 720 bp

Locus ID: 21949


UniProt ID: [P32972](#)

Cytogenetics: 4 34.06 cM

MW: 26.5 kDa

Gene Summary: Cytokine that binds to TNFRSF8/CD30. Induces proliferation of T-cells.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219380