

Product datasheet for **MR219382L4V**

Mapre2 (NM_001162941) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mapre2 (NM_001162941) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mapre2
Synonyms:	AI314113; C820009F03Rik; D18Abb1e; EB1; EB2; RP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001162941
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219382).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162941.1 , NP_001156413.1
RefSeq Size:	3993 bp
RefSeq ORF:	954 bp
Locus ID:	212307
UniProt ID:	Q8R001
Cytogenetics:	18 12.08 cM
Gene Summary:	May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)