

Product datasheet for **MR219394**

Gimap4 (NM_175048) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gimap4 (NM_175048) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gimap4
Synonyms: AU019574; E430007K16Rik; I; lan-1; lan1; lmap4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR219394 representing NM_175048
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAAGTCCAGTGCGGTGGTGCGGGTT**CATCCAGAAAGTTCAAGGAGCAGCCATGAGCTT**GAAACC
AAGATCAAGGAAT**CCCCAACTGAGAATTGTCTTACTTGAAAACTGGAGCAGGAAAGAGTTCAACAGG**
GAACAAGGTGATGCATGAGTTCCAGAATCGCTACTGTTTGTTCACAACAGAGCCTCAGGCGCTGAAAAG
GAAGAGCAGAAGATGCAATTGTTGACCTTGGTCCAGAGCATGGATGAAGTAATGAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219394 representing NM_175048
Red=Cloning site Green=Tags(s)

MEVQCGGAGFIPSSRSSELGNQDQGIPLRIVLLGKTGAGKSSTGNKVMHEFQNYRCLFNNRASGAEK
EEQKMQLLTLVQSMDEVK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_175048

ORF Size: 267 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175048.3</u> , <u>NP_778213.2</u>
RefSeq Size:	967 bp
RefSeq ORF:	270 bp
Locus ID:	107526
Cytogenetics:	6 B2.3
MW:	10.3 kDa
Gene Summary:	This gene encodes a protein belonging to the GTP-binding superfamily and to the immun-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. This gene exists within a cluster of other related genes located on mouse chromosome 6. This family member encodes a lymphoid signaling protein that functions to accelerate programmed T-cell death, which appears to correlate with the phosphorylation status of the protein. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]