

Product datasheet for **MR219403A1V**

Mouse Il36rn (NM_019451) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Il36rn (NM_019451) AAV Particle
Tag: Myc-DDK
Symbol: Il36rn
Synonyms: AI413231; Fil1delta; Il-1h3; IL-36Ra; Il1hy1; IL36RN
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR219403 representing NM_019451
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGTTCTGAGTGGGGCACTATGCTTCCGAATGAAGGATTCAGCCTTGAAGGTACTGTATCTGCACAATA
 ACCAGCTGCTGGCTGGAGGACTGCACGCAGAGAAGGTCATTAAAGGTGAGGAGATCAGTGTGTCCCAAA
 TCGGGCACTGGATGCCAGTCTGTCCCTGTCTCCTGGGCGTTCAAGGAGGAAGCCAGTGCCTATCTTGT
 GGGACAGAGAAAGGGCCAATTCTGAACTTGAGCCAGTGAACATCATGGAGCTCTACCTCGGGGCCAAGG
 AATCAAAGAGCTTCACCTTCTACCGCGGGATATGGGTCTTACCTCCAGCTTCGAATCCGCTGCCTACCC
 AGGCTGGTTCTCTGCACCTCACCGGAAGCTGACCAGCCTGTCAGGCTCACTCAGATCCCTGAGGACCCC
 GCCTGGGATGCTCCCATCACAGACTTCTACTTTCAGCAGTGTGAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219403 representing NM_019451
 Red=Cloning site Green=Tags(s)

MVLSGALCFRMKDSALKVLYLHNNQLLAGGLHAEKVIKGEESVVPNRALDASLSPVILGVQGSQCLSC
 GTEKGPIKLPEVNIMELYLGAKESKSFYRRDMGLTSSFESAAYPGWFLCTSPEDQPVRLTQIPEDP
 AWDAPITDFYFQQCD

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2


[View online »](#)

ACCN:	NM_019451
ORF Size:	465 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_019451.2</u>
RefSeq Size:	1759 bp
RefSeq ORF:	471 bp
Locus ID:	54450
UniProt ID:	<u>Q9QYY1</u>
Cytogenetics:	2 16.31 cM
MW:	17 kDa