

## Product datasheet for **MR219403L4V**

### II1f5 (NM\_019451) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	II1f5 (NM_019451) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	II1f5
Synonyms:	AI413231; Fil1 delta; Il-1 h3; IL-36Ra; Il1hy1; IL36RN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019451
ORF Size:	468 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219403).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019451.2</a>
RefSeq Size:	1759 bp
RefSeq ORF:	471 bp
Locus ID:	54450
UniProt ID:	<a href="#">Q9QYY1</a>
Cytogenetics:	2 16.31 cM



[View online »](#)

**Gene Summary:**

Inhibits the activity of interleukin-36 (IL36A,IL36B and IL36G) by binding to receptor IL1RL2/IL-36R and preventing its association with the coreceptor IL1RAP for signaling. Part of the IL-36 signaling system that is thought to be present in epithelial barriers and to take part in local inflammatory response; similar to the IL-1 system with which it shares the coreceptor. Proposed to play a role in skin inflammation. May be involved in the innate immune response to fungal pathogens. May activate an anti-inflammatory signaling pathway by recruiting SIGIRR.[UniProtKB/Swiss-Prot Function]