

Product datasheet for **MR219429A1V**

Mouse Il36g (NM_153511) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Il36g (NM_153511) AAV Particle
Tag: Myc-DDK
Symbol: Il36g
Synonyms: Il36g
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR219429 representing NM_153511
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTTTTCTAAACACCCATTTTCTACACACATCTCAGGAAGAGAACTCCTGACTTTGGGGAGGTTTTTG
ACTTGGACCAGCAGGTGTGGATCTTTCGTAATCAGGCCCTTGACAGTTCCACGAAGCCACAGAGTAAC
CCCAGTCAGCGTGAATCCTCCCATGCAAGTACCCAGAGTCTTGAACAGGACAAAGGGATTGCCATT
TATTTGGGAATTCAGAATCCAGATAAATGCCTGTTTTGTAAGGAAGTTAATGGACACCCTACTTTGCTGC
TAAAGGAAGAGAAGATTTGGATTTGTACCACCACCTGAGCCAATGAAGCCATTCTGTTTTACCACAC
CCGGACAGGTGGAACATCCACCTTTGAATCAGTGGCTTTCCCTGGCCACTATATTGCCTCCTCCAAGACT
GGCAACCCCATCTTCCTCACATCAAAAAAGGGAGAATATTACAACATTAATTCAATTTAGATATAAAGT
CT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219429 representing NM_153511
Red=Cloning site Green=Tags(s)

MFSKHPFSTHISGRETPDFGEVFDLDQQVWIFRNQALVTVPRSHRVTPVSVTILPCKYPESLEQDKGIAI
YLGIQNPDKCLFCKEVNGHPTLLLKEEKILDLYHHPEPMKPFIFYHTRTGGTSTFESVAFPGHYIASSKT
GNPIFLTSSKGEYNNINFNLDIKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2



[View online »](#)

ACCN:	NM_153511
ORF Size:	492 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_153511.3</u>
RefSeq Size:	1647 bp
RefSeq ORF:	582 bp
Locus ID:	215257
Cytogenetics:	2 16.24 cM
MW:	22.5 kDa